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CE  FC

PB

Designed in Hong Kong SAR
www.zikkotech.com

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Dr. Rock Mini

BianStone far infrared heating therapy massager

Product Introduction

This revolutionary massager adopts a brand-new heat energy injection technology, which will give you an incomparable feeling of total soothing and relaxation from the first use. It can heat up extremely quickly, emit far infrared rays into the body, and produce the same frequency resonance with cells, activate cells, promote blood circulation, relieve pain, relax muscles, and restore physical fitness. The heat energy supports the vertical skin-penetrating ultrasonic wave, which quickly injects the effective essential oils into the skin below the epidermis, enhances the skin's brilliance, and makes the skin smoother. You can also download the App through your mobile phone to create and play your own unique physiotherapy program. Enjoy the smart warm healthy care of pure natural BianStone (Scraping Stone) in your palm.

Product schematic

control panel
USB-C charging port
Physiotherapy head (5 replaceable physiotherapy heads)
Heating key
Power ON/OFF
Vibrate key
Charging indicator
Vibration indicator
Heating indicator

Installation method and characteristics

- The installation of the massage head is straight plug, no need to twist.
- The red line of the Bianstone physiotherapy head is aligned with the piston port.
- Need to make sure to insert to the bottom.

Bian Stone semicircle heating therapy head
40~60 °C temperature
Far-infrared physiotherapy relieves sciatica, frozen shoulder, joint pain, muscle aches, swelling, meridian dredging and other pains, lymphatic drainage, eliminates dark circles and eye bags, warms the belly (regulates the stomach and intestines), and quickly introduces essential oils (beauty skin) Body).

Ball-shaped head
It is suitable for massaging large muscle groups such as arms, waist, back, buttocks, thigh, calf, and etc.

U-shaped head
It is suitable for massaging neck, spine and achilles tendon.

Cylindrical head
Suitable for impacting deep tissues, such as muscles, palms, soles, etc.

Flat head
Suitable for relaxation and plasticity of various muscle parts.

Charging instructions

- Before the first use, charge it for about 4 hours and fully charge it;
- When charging, a dedicated USB cable will be provided to connect to the USB-C socket of the device, and any 5V charger head will be used on the other end. (5V/2A recommended);
- All power indicators on the instrument flash regularly, indicating that it is charging;
- When all the battery indicator lights on the instrument stop flashing and stay on, it means that the battery is fully charged;
- The battery can be charged at any time and at any power level under specified environmental conditions;
- When the power is less than 25%, it should be charged in time;
- It is not recommended to completely drain the battery to 25%;

The charging port is located at the top of the product handle

8. The average running time of the equipment after loading is three hours, which depends on the vibration frequency and applied pressure during use.

Instructions for use

Heating key Vibration key

- Turn on, vibrate:**
Long press the button for 2 seconds, when the vibration is turned on, it will be level 1, and the display light will light up. Short press once for 2 levels, and two display lights will light up. Short press again for 3 levels, and the display lights will turn on three. Continue to short press to turn off the vibration function, and all three display lights will go out;
- Heating:**
Press and hold the button for 2 seconds, the hot compress is turned on, the temperature is about 40°C, and the heating display light turns on. Short press once, the temperature is about 50°C, and two display lights will light up. Short press again, the temperature is about 60°C, and three display lights will light up. Continue to short press once to turn off the hot compress function, and all three display lights will go out.
- Note:**
When either function of vibration or hot compress is turned on, short press to turn on another function, vibration gear or heating temperature adjustment, refer to 1 or 2;
- Bianstone therapy only**
With 2 hot compress functions, other physiotherapy heads have no hot compress function when turned on;
- Independent use:**
Both vibration and hot compress functions can be used independently;
- 15 minutes automatic shutdown:**
The automatic timing of vibration and hot compress function is 15 minutes, and it will automatically turn off after 15 minutes of use;
- Manual shutdown:**
Long press any key in the power-on state to turn off the power.

Dr.Rock APP operating instructions

- The Dr.Rock App is available on the Apple App Store. Please search for keywords such as , "Dr.Rock", etc. on the App Store to download;

2. After downloading and installing, turn on the phone's Bluetooth → click "Dr.Rock APP" → "Device" → "Scan Bluetooth", the product ID will pop up, and click Connect → "Bind";

3. After the product is connected, adjust and use it reasonably according to your own needs and the option functions on the APP;

Vibration speed Celsius display
Bluetooth switch Bind/unbind switch
Default automatic mode Custom mode switch (manual adjustment)
Vibration adjustment Heating switch
Saved record Celsius adjustment

4. **Introduction to other functions of Dr. Rock APP**
① "Products" → You can browse and understand other Dr. Rock series products, and you can also purchase products on the APP;
② "Physiotherapy" → Learn how to use Dr.Rock series products.

Maintenance, cleaning, storage

- After using the product, clean it with a damp towel, and then dry it with a soft cloth;
- Before and after using the Bianstone physiotherapy head, wipe and disinfect the surface of the Bianstone with medical alcohol;
- Turn off the power after each use of the product;
- When the product is stored for a long time or travels, please put it in the special storage bag provided;
- When the product is stored for a long time, the product should be charged once a month to avoid damage to the battery or moisture;
- Pay attention to waterproof, fireproof, damp and drop when placing the product. At the same time, prevent minors from using it.
- If the product is used for a long time, in order to maintain good noise reduction, it is necessary to add lubricant between the sleeve ports of the machine body (as shown in the right figure)

Lubricant
muzzle

Specifications

Model : H-MG400
Name : Dr.Rock mini BianStone far infrared heating therapy massager
Impact stroke : 6mm / Strength : 8KG
Strength gear : 3 gears variable speed adjustment
Gear/speed : 1 gear 1000 rpm, 2 gear 2000 rpm, 3 gear 3000 rpm
Hot compress temperature : 1 gear 40°C, 2 gear 50°C, 3 gear 60°C
Charging input : USB-C 2A/5V
Battery : 7.4V/2200mAh, 16.28Wh, lithium battery
Endurance time : 2-4 hours (varies according to usage status)
Charging time : about 4 hours
Size : 164*48*158mm
Weight : 492g
Color : gray

Packaging accessories

- Dr.Rock Mini x1
- User's manual x1
- Charging cable x1
- Storage bag x1
- Ball massage head x1
- Bian Stone semicircle head x1
- Flat massage head x1
- U shape massage head x1
- Cylindrical massage head x1

Bian Stone far infrared warm compress physiotherapy

Warm compress therapy is an important method in stone therapy. Its origin is to use baked hot stone to rest on the affected part of the disease or massage, which can relieve pain and treat the disease. This therapy is called stone warm compress. Electrothermal stone therapy is a kind of warm compress, i.e. current is used to heat the stone through a special instrument. The far infrared electromagnetic wave radiation heat energy generated by the heated stone can be conducted in channels and collaterals to achieve the purpose of treating various diseases. It is a clever combination of traditional stone therapy and modern high technology. At the same time, stone has a strange energy field, which acts on the surface of human skin to generate extremely far infrared radiation with extremely wide frequency band and far infrared frequency ranging from 7 to 20 microns. Every time stone rubs the human body, it can generate ultrasonic pulses beneficial to health. The average number of ultrasonic pulses can reach 3706, with the frequency range of 2-2 million Hz. Therefore, stone is the first of all materials.

Dr.Rock series Bian Stone far infrared warm compress physiotherapy has similar effect to that of warm compress therapy. It is a compound therapy of stone hyperthermia and massage, and has excellent curative effects of promoting blood circulation, removing blood stasis and dredging channels and collaterals. After thousands of years of practice, warm compress therapy has proved to have a wide range of health care functions. Clinical evidence shows that electrothermal stone therapy has good effects of warming channels and dredging veins, eliminating arthralgia, relieving spasm and sedation, and can be used for various pains, relieving muscle spasm, digestive system diseases, and etc. It also shows great advantages in beauty treatment, weight loss, insomnia adjustment and fatigue elimination.

Massage impact therapy

Impact on skin and subcutaneous tissue

Loosening thickened connective tissue and fascia, breaking the barrier of releasing fascia contraction, removing adhesion and internal scar tissue after injury or surgery, increasing venous return and lymphatic circulation, and improving muscle extensibility. Especially when muscle spasm, stiffness and restriction occur, the reflex contraction of muscle fibers is punctured, and the range of motion is increased by relieving stiffness of joints.

Influence on central and sympathetic nervous system

Stimulate local muscle receptors to relieve pain and help relax, without causing pain or inflammation.

Influence on muscle, preserved muscle, arthralgia

Common causes of muscle pain are trauma, overuse, and tension. After injury or trauma, impact therapy can speed up the healing of the inflamed area. Similarly, when a part of the body is overused or stressed, impact therapy can speed up blood circulation in the damaged part and reduce potential tension.

Exercise rehabilitation

Physiotherapy helps the body recover from strenuous exercise, and is beneficial

for conditioning and rebuilding training. Rehabilitation will remove metabolic waste and better nourish cells, normalize tissue elasticity and speed up wound healing.

Recovery after operation and injury

After trauma or surgical repair, the stimulation of impact therapy can eliminate faults and avoid formation and adhesion of scar tissue during the repair process of inflammatory activity of the body.

Relax

Impact physiotherapy can repair the body, increase the patient's body comfort and improve the overall well-being

indications and contraindications

Dr.Rock series deep penetration can relieve pain, bursitis, muscle spasm, neuralgia and other musculoskeletal diseases caused by joint sprain in a short time. Almost all soft tissue diseases are related to trauma, overuse of muscles or abnormal posture. Dr.Rock can be used to treat scar tissue involved in most of these diseases. When using Dr.Rock, excessive pressure and long-term physical therapy should be avoided. Sensitive areas such as head, face and cervical vertebra should also be avoided. Physical therapy should be stopped immediately if there is inflammation, swelling or pain aggravation. The current research is about the impact of vibration training on specific medical diseases. This is likely to reduce the types of contraindications mentioned above. Practical experience shows that in many cases, the combination of vibration training with physiotherapy program is also a physiotherapy suggestion. This must be done on the advice or accompanied by a doctor, expert or physiotherapist.

Indication

- ✓ Pain and spasm due to muscle contusion, sprain and strain
- ✓ It can help edema fluid flow in swollen parts,
- ✓ Relax thickened connective tissue and fascia,
- ✓ Reduce lactic acid accumulation in muscles, increase joint activity,
- ✓ Eliminate muscle fatigue and
- ✓ Sciatica (Bian stone hot compress physiotherapy) and
- ✓ Scapulohumeral periarthritis (Bian stone hot compress physiotherapy)
- ✓ Periarthritis (Bian stone hot compress physiotherapy)

Contraindications

- Aneurysm, hemorrhage, heart disease with blood thinner;
- cardiac pacemaker or defibrillator;
- Pregnancy, cancer;
- It cannot be used within 3 inches of internal fixation when joint internal fixation less than 90 days,
- Sensitive areas: head, face, cervical spine, Spinal vertebrae
- Muscles cover less locations near bones:
Tibia, back of foot, back of hand

Physiotherapy considerations

After injury

- The first physical therapy after injury should be based on the patient's satisfaction and comfort.
- ROM evaluation
- Determine the location and surrounding muscles according to the physiotherapy plan
- Start the instrument and set it at low to medium intensity
- The pressure applied is the weight of the instrument itself

After surgery

- The first postoperative physical therapy gradually increases the time and intensity of physical therapy based on the patient's response and feelings
- Determine the treatment site
- Ask the patient if they have open wounds, stitches, drainage, and internal fixation. If yes, do not irritate the area
- Start the instrument and set it at the lowest level
- Cover the treatment area with a towel
- Use a physiotherapy device to apply pressure on the physiotherapy area
- Physiotherapy is limited to the area around the surgical site

Lymphatic drainage

- Requires mild, systemic physical therapy
- Slowly move from the distal end to the proximal end of the limb during physical therapy, that is, in the direction of the return of the heart's veins

Effect

- Reduce swelling and pain
- Eliminate scar tissue after surgery and injury
- Restore muscle elasticity and extensibility
- Trigger point therapy
- Promote blood circulation and lymphatic return
- Myofascial relaxation

Note : Whether it is controlled depression or applying pressure, bruises may occur. Check the physiotherapy area frequently, and stop using it as soon as there is pain or discomfort.

Safety Precautions

- It is for adult use only, and it is forbidden for minors to use it;
- It can only be used on the dry and clean skin surface of the body part. It can also be used through clothes. Gently press and move on the skin for about 60 seconds for each part;
- It is forbidden to use fatal acupuncture points on the head or body that cannot bear any force, and any hard and boned parts;
- It is forbidden to immerse the product in water or put it in a fire, close to the fire source or high temperature;

- It is forbidden to drop the product or use it improperly;
- It is forbidden to touch the shaft and back of the vibrating head with fingers, hair or other parts of the body to avoid being caught;
- It is forbidden to change or modify the product in any way;
- It is forbidden to let the product run or charge unattended;
- Before each use, please carefully check whether all the functions of the product are in normal use;
- The product Bianstone physiotherapy head has a hot compress function. Please adjust it according to your own needs before use to avoid burns. People who are not sensitive to heat should use it with caution.

Abnormal faults and troubleshooting

Fault description	Abnormal reason	Exception handling measures
The product cannot be turned on	Check whether the battery is sufficient	Please charge (please charge for more than 6 hours if not used for a long time)
Automatic shutdown during product use	1.Low battery 2.The product has been working continuously for 15 minutes 3.The product is protected by excessive resistance during use	1. Please charge 2. Long-term protection mode, just restart 3. Normal protection, just restart.
Product cannot be charged	1. Check whether the port is well connected and whether the adapter supports 5V/2A 2. The product temperature or ambient temperature above 40°C cannot be charged.	1. Please re-plug or replace the matching adapter 2. Please wait for the product to cool down or recharge in a normal temperature environment.
Physiotherapy head slipped and dropped	Check whether the physiotherapy head is installed in place	Operation reference manual physiotherapy head installation method
Bianstone head temperature is too high	The hot compress mode is turned on in 3 levels. If you have not touched the body for a long time, it will balance the temperature of the Bianstone	Turn off the hot compress mode and wait for the temperature of the Bianstone to cool down or lower the hot compress mode to use
Excessive noise during product use	The partial silencer grease of the product is dry or caused by other reasons	For specific reasons, please contact after-sales processing

Note : If there are other abnormalities, please contact customer service

After-sales service

One year warranty

The whole machine is guaranteed for one year from the date of purchase, and the attached accessories are guaranteed for 6 months. If the quality problems caused by product defects cannot be used normally, our company will be responsible for free maintenance, replacement, or replacement of necessary components.

Note: The following disclaimer is not within the scope of this warranty.

- Damage caused by accident, abuse, improper handling or transportation;
- Repair the product without permission;
- The product is not maintained or used according to the instructions;
- The loss is higher than the cost of the product;
- Quality problems in delivery due to abnormal storage or protection conditions at the customer's workplace;
- Unable to provide valid proof of purchase;
- Products purchased from unofficial or authorized channels

Some states and countries do not allow restrictions on fault compensation, so the above restrictions may not be available to you. This warranty covers specific legal rights, but other rights may vary by country and state. This warranty service is provided by Zikko (Hong Kong) Digital Technology Co., Ltd. and Zikko (China) Co., Ltd. If you are a customer, please contact the physical store, distributor or retailer where you purchased the product directly to learn about the warranty and return procedures.

Environmentally friendly recycling

The product has a built-in rechargeable battery, which may contain substances that pollute the environment. Before discarding the product, be sure to turn off the power and take out the battery, or send the product to an officially designated battery recycling place. You can also send the discarded products back to our designated recycling department (please contact after-sales for recycling information), and the maintenance department will decompose and recycle the battery.

- Do not short-circuit the battery to prevent the battery from burning and causing a fire.
- Do not throw the battery into the fire to prevent explosion or release of toxic substances

Name and content of hazardous substances in this product

Part Name	Toxic and harmful substances or elements					
	Lead (Pb)	Mercury (Hg)	Cadmium (Cd)	Hexavalent chromium (Cr(VI))	PBBs (PBB)	PBDEs (PBDE)
plastic shell	O	O	O	O	O	O
metallic parts	O	O	O	O	O	O
Rubber parts	O	O	O	O	O	O
Circuit board	O	O	O	O	O	O
Connection line	O	O	O	O	O	O
Motor	O	O	O	O	O	O

This form is compiled in accordance with SJ/T 11364
O means that the content of the hazardous substance in all homogeneous materials of the part is below the limit requirement specified in GB/T 26572.
X indicates that the content of the hazardous substance in at least one of the homogeneous materials of the part exceeds the limit requirement specified in the GB/T26572 standard

1. During the environmental protection period, please use this product in accordance with the instruction manual.
2. This environmental protection use period does not cover vulnerable parts: batteries.

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