

Dr.Rock

Energy Pain Relief Patch

能量鎮痛貼



Scan QR Code
請掃描QR碼



Dr.Rock
Designed in Hong Kong SAR
www.drrock.com

User Manual

使用說明書

After Service Card

售後服務保障卡

ENGLISH 1 ~ 10
繁體中文 11 ~ 19

Dr.Rock 中文頻道
DrRock English channel




掃碼觀看視頻教程
Scan QR Watch usage video tutorial

Wish you health and peace.
願健康與平安歸給您。



Product Name
產品名稱

Product Model
產品型號

Product Serial No.
產品序號

Purchasing Date
購買日期

Customer Name
顧客姓名

Phone
聯系電話

E-Mail
郵箱

Addr.
聯系地址

Distributor Name
經銷商名稱

Phone
聯系電話

E-Mail
郵箱

Sign/Seal Here:
經銷商簽章

Dr.Rock Energy Pain Relief Patch

Product Description

Dr.Rock Energy Pain Relief Patch is a non-invasive therapeutic device that can be magnetically attached and detached from the main unit. It utilizes graphene heating technology to release far-infrared rays and negative ions. The patch also contains PEMF (Pulsed Electromagnetic Field) energy at 7.83Hz, along with a natural trace element formula in the energy patch, to achieve pain relief effects.

Six Key Features

- Graphene heating
- PEMF 7.83Hz
- 30-45°C thermostatic
- Vibration adjustment
- Magnetism
- Easy to carry

Product Diagram



Product User Manual

USB-C Port
Temp Indicator Light
Temp Knob
Power BTN / Vibration BTN / PEMF BTN
Battery indicator light
heating indicator light
Vibration mode indicator light
LED Screen

Loc & Func

Power BTN / Vibration BTN / PEMF BTN
Temp Knob
Temp Indicator Light
Vibration mode indicator light
heating indicator light
Battery indicator light

Instructions for use

Long press for 2 seconds to power on/off. After powering on, short press to adjust vibration intensity. Adjustment sequence: (1st gear) (2nd gear) (3rd gear) - off. After powering on, PEMF (Pulsed Electromagnetic Field) is turned on by default, with a frequency of 7.83Hz. Quick double-click to enable/disable PEMF function. Use the temperature control knob to adjust the temperature within the range of 30 to 45°C. Clockwise rotation increases the temperature by 1 degree until it reaches 45°C, then slope. Counterclockwise rotation decreases the temperature by 1 degree until it reaches 30°C, then stops. When selecting the minimum temporary, counterclockwise rotation turns off the heat therapy. When the heat therapy is turned off, rotating clockwise turns it on. "30-45°C" indicates the current heat therapy temperature setting. "00" means heat therapy is turned off. The flashing of the "i" icon indicates that heating is in progress. When the set temperature is reached, the "i" icon disappears, indicating that the heating operation has stopped. Display the current battery power status. Full battery: 100% to 80% - Battery indicator light is green. Normal battery level: 80% to 30% - Battery indicator light is blue. Low battery level: 30% to 0% - Battery indicator light is red. When the battery is low, there will be a "beep" sound notification. Charging status: Red breathing effect on the battery indicator light. Fully charged: Battery indicator light remains solid green. This product has a memory function, which retains the previous temperature setting even after power off or separation from the controller.

Product Usage Diagram



1. Open the packaging of the energy patch and remove the patch.
2. Peel off the entire adhesive film with the Dr.Rock logo.
3. Stick the patch onto the device.
4. Peel off the film and keep it until it runs out, then stick it back on for the next use.
5. Long press the power button to turn on the device.
6. Gently press and apply it to the skin, adjusting the placement according to different areas of pain.
7. When adjusting the temperature or function, you can detach the controller separately for adjustment.
8. Place the controller gently to magnetically attach it to the pain relief patch.
9. After use, remove the controller first.
10. Then remove the pain relief patch from the skin.
11. After turning off the device, place it back in the product storage pouch for future use.
12. After turning off the device, place it back in the product storage pouch for future use.

1. Use of Joint pain energy patch

1. Open the Energy Patch package and take out the patch. (Color is ●)

2. Remove all of the anti-adhesive film and apply the patch to the heating film.

3. Stick it on the skin and use it according to the different parts of the controller on the Energy Pain Relief Patch.

2. Use of Waist and leg pain energy patch

1. Open the Energy Patch package and take out the patch. (Color is ●)

2. Remove all of the anti-adhesive film and apply the patch to the heating film.

3. Stick it on your skin and use it according to the different parts of your neck and shoulder pain.

4. Open the controller to set the desired temperature and gently place the controller on the Energy Pain Relief Patch.

3. Use of Shoulder neck pain energy patch

1. Open the Energy Patch package and take out the patch. (Color is ●)

2. Remove all of the anti-adhesive film and apply the patch to the heating film.

3. Stick it on your skin and use it according to the desired temperature and gently place the controller on the Energy Pain Relief Patch.

4. After use, remove the controller first.

5. Put on the Protective sleeve to prevent it from falling off.

6. After use, remove the Protective sleeve.

4. Use of Neck and shoulder pain energy patch

1. Open the Energy Patch package and take out the patch. (Color is ●)

2. Remove all of the anti-adhesive film and apply the patch to the heating film.

3. Stick it on your skin and use it according to the desired temperature and gently place the controller on the Energy Pain Relief Patch.

4. After use, remove the controller first.

5. Put on the Protective sleeve to prevent it from falling off.

6. After use, remove the Protective sleeve.

5. Use of Neck and shoulder pain energy patch

1. Open the Energy Patch package and take out the patch. (Color is ●)

2. Remove all of the anti-adhesive film and apply the patch to the heating film.

3. Stick it on your skin and use it according to the desired temperature and gently place the controller on the Energy Pain Relief Patch.

4. After use, remove the controller first.

5. Put on the Protective sleeve to prevent it from falling off.

6. After use, remove the Protective sleeve.

Specifications

Model: EP100
Name: Dr.Rock Energy Pain Relief Patch
Auxiliary heating: graphene membrane
Hot compress temperature: 30-45°C
Adjustable temperature: 30-45°C (Note: The above will vary by ±5°C due to environment and usage)
Vibration: 3-speed adjustment
Automatic timing: every 20 minutes. (To use again, please remove the controller and restart.)
Ultra long wave frequency: 7.83 Hz
Power input: USB-C 5V/1.2A
Lithium battery: 3.7V/1500mAh, 5.5Wh
Battery life: about 1 hour (depending on usage status)
Charging time: about 2 hours
Size: 149*79*30mm / Weight: 88g / Colour: Black

FAQs & Troubleshooting

Fault description	Reason for the anomaly	Handling measures
The product cannot power on	Check if the battery is fully charged.	Please charge the device. If it has not been used for a long time, please charge it for at least 2 hours.
The device shuts down during use	Automatic timing has ended or insufficient battery power.	Please reboot or charge again.
The product cannot charge or shows no charging indication	1. Check if the charging port is properly connected. 2. Check if the charging power source is sufficient power.	1. Please reinsert or release the appropriate adapter. 2. Check if the charging power source is working properly.
The product doesn't vibrate	The vibration motor is not receiving power or is damaged.	For specific issues, please contact customer support for further assistance.
The product doesn't feel warm	Wait for 1-2 minutes. The heating element may not be receiving power or it could be burned out.	For specific issues, please contact customer support for further assistance.

Note: If there are any other abnormal conditions, please contact customer service promptly.

After-Sales Service

One-year warranty

From the date of purchase, the whole machine is covered by a one-year warranty, while the accompanying accessories are covered by a six-month warranty. In the event of quality issues caused by manufacturing defects that prevent normal usage, our company will be responsible for free repairs, replacements, or necessary component replacements.

Note: The following disclaimer is not applicable within this warranty scope:

- Damages caused by accidents, misuse, mishandling, or transportation;
- Unauthorized repairs of the product;
- Failure to follow the maintenance or usage instructions stated in the user manual;
- Loss exceeding the cost of the product;
- Quality issues arising from abnormal storage or protection conditions at the customer's workplace;
- Inability to provide a valid purchase receipt;
- Products purchased from unofficial or unauthorized channels;

Environmental Protection And Recycling

The product contains a built-in rechargeable battery that may contain environmentally hazardous substances. Before disposing of the product, please make sure to turn off the power and remove the battery, or take the product to an official designated battery recycling center. You can also send the discarded product to our designated recycling department (please contact customer service for recycling information), where the battery will be disassembled and recycled by the repair department.

Please do not short-circuit the battery to prevent battery combustion and fire hazards.

Do not dispose of the battery in fire to prevent explosions or the release of toxic substances.

Safety Precautions

1. Before use, carefully read the instruction manual and follow the operating instructions.
2. Do not use the product alone by individuals who are unable to take care of themselves, infants, or those who are not sensitive to heat.
3. When using, pay attention to the temperature after heating and select the appropriate temperature setting that is tolerable for your skin. Avoid overheating to prevent burns.
4. It is strictly prohibited to use chemical solvents to clean the device body, charging cable, etc.
5. Do not dismantle the device by yourself in any way, otherwise, any resulting damage will be solely your responsibility.
6. If you feel the temperature is too high or if the applied area becomes red, it is recommended to lower the heat setting or temporarily discontinue use.

Maintenance & Storage

1. After using the product, remove the used energy patch promptly.
2. When storing the product, keep it in a dry, well-ventilated place, protected from moisture, fire, sunlight, pressure, impact, and contamination.
3. Once the energy patch is removed from its packaging, it should be used immediately and should not be placed back for storage to avoid contamination and loss of adhesive properties.
4. After using the main unit of the product, turn off the power and place it in a safe location, away from children's reach.
5. When storing the product for a long time, please expose the product to direct sunlight, as it may cause fading, deformation, or even burning. Additionally, if not used for a long time, it is recommended to charge the main unit of the product once a month to prevent battery damage or moisture issues.
6. Do not use this product for any other purposes.

Dr.Rock 能量鎮痛貼

產品簡介

Dr.Rock能量鎮痛貼是一款主機可磁吸分離、非侵入式的理療設備。通過石墨烯加熱，釋放遠紅外線和負離子，內含PEMF（脈衝磁場）7.83Hz能量，配合使用能量貼的天然微量元素配方，達到舒緩疼痛效果。

六大特點

- 石墨烯加熱
- 脈衝磁場7.83Hz
- 30-45°C電溫
- 震動調節
- 分離式貼敷
- 輕巧易帶

產品示意圖



產品使用說明

USB-C充電口
溫度設置指示燈
調整按鈕
PEMF指示燈
電源鍵/震動鍵/PEMF鍵
電池顯示燈
加熱指示燈
震動指示燈
LED顯示屏

鍵位和功能

長按2s為開/關機條件。開機後，短按為開啟震動強度調節：調節順序：(檔1) (檔2) (檔3) (檔4) (檔5) (檔6) (檔7) (檔8) (檔9) (檔10) (檔11) (檔12) (檔13) (檔14) (檔15) (檔16) (檔17) (檔18) (檔19) (檔20) (檔21) (檔22) (檔23) (檔24) (檔25) (檔26) (檔27) (檔28) (檔29) (檔30) (檔31) (檔32) (檔33) (檔34) (檔35) (檔36) (檔37) (檔38) (檔39) (檔40) (檔41) (檔42) (檔43) (檔44) (檔45) (檔46) (檔47) (檔48) (檔49) (檔50) (檔51) (檔52) (檔53) (檔54) (檔55) (檔56) (檔57) (檔58) (檔59) (檔60) (檔61) (檔62) (檔63) (檔64) (檔65) (檔66) (檔67) (檔68) (檔69) (檔70) (檔71) (檔72) (檔73) (檔74) (檔75) (檔76) (檔77) (檔78) (檔79) (檔80) (檔81) (檔82) (檔83) (檔84) (檔85) (檔86) (檔87) (檔88) (檔89) (檔90) (檔91) (檔92) (檔93) (檔94) (檔95) (檔96) (檔97) (檔98) (檔99) (檔100) (檔101) (檔102) (檔103) (檔104) (檔105) (檔106) (檔107) (檔108) (檔109) (檔110) (檔111) (檔112) (檔113) (檔114) (檔115) (檔116) (檔117) (檔118) (檔119) (檔120) (檔121) (檔122) (檔123) (檔124) (檔125) (檔126) (檔127) (檔128) (檔129) (檔130) (檔131) (檔132) (檔133) (檔134) (檔135) (檔136) (檔137) (檔138) (檔139) (檔140) (檔141) (檔142) (檔143) (檔144) (檔145) (檔146) (檔147) (檔148) (檔149) (檔150) (檔151) (檔152) (檔153) (檔154) (檔155) (檔156) (檔157) (檔158) (檔159) (檔160) (檔161) (檔162) (檔163) (檔164) (檔165) (檔166) (檔167) (檔168) (檔169) (檔170) (檔171) (檔172) (檔173) (檔174) (檔175) (檔176) (檔177) (檔178) (檔179) (檔180) (檔181) (檔182) (檔183) (檔184) (檔185) (檔186) (檔187) (檔188) (檔189) (檔190) (檔191) (檔192) (檔193) (檔194) (檔195) (檔196) (檔197) (檔198) (檔199) (檔200) (檔201) (檔202) (檔203) (檔204) (檔205) (檔206) (檔207) (檔208) (檔209) (檔210) (檔211) (檔212) (檔213) (檔214) (檔215) (檔216) (檔217) (檔218) (檔219) (檔220) (檔221) (檔222) (檔223) (檔224) (檔225) (檔226) (檔227) (檔228) (檔229) (檔230) (檔231) (檔232) (檔233) (檔234) (檔235) (檔236) (檔237) (檔238) (檔239) (檔240) (檔241) (檔242) (檔243) (檔244) (檔245) (檔246) (檔247) (檔248) (檔249) (檔250) (檔251) (檔252) (檔253) (檔254) (檔255) (檔256) (檔257) (檔258) (檔259) (檔260) (檔261) (檔262) (檔263) (檔264) (檔265) (檔266) (檔267) (檔268) (檔269) (檔270) (檔271) (檔272) (檔273) (檔274) (檔275) (檔276) (檔277) (檔278) (檔279) (檔280) (檔281) (檔282) (檔283) (檔284) (檔285) (檔286) (檔287) (檔288) (檔289) (檔290) (檔291) (檔292) (檔293) (檔294) (檔295) (檔296) (檔297) (檔298) (檔299) (檔300) (檔301) (檔302) (檔303) (檔304) (檔305) (檔306) (檔307) (檔308) (檔309) (檔310) (檔311) (檔312) (檔313) (檔314) (檔315) (檔316) (檔317) (檔318) (檔319) (檔320) (檔321) (檔322) (檔323) (檔324) (檔325) (檔326) (檔327) (檔328) (檔329) (檔330) (檔331) (檔332) (檔333) (檔334) (檔335) (檔336) (檔337) (檔338) (檔339) (檔340) (檔341) (檔342) (檔343) (檔344) (檔345) (檔346) (檔347) (檔348) (檔349) (檔350) (檔351) (檔352) (檔353) (檔354) (檔355) (檔356) (檔357) (檔358) (檔359) (檔360) (檔361) (檔362) (檔363) (檔364) (檔365) (檔366) (檔367) (檔368) (檔369) (檔370) (檔371) (檔372) (檔373) (檔374) (檔375) (檔376) (檔377) (檔378) (檔379) (檔380) (檔381) (檔382) (檔383) (檔384) (檔385) (檔386) (檔387) (檔388) (檔389) (檔390) (檔391) (檔392) (檔393) (檔394) (檔395) (檔396) (檔397) (檔398) (檔399) (檔400) (檔401) (檔402) (檔403) (檔404) (檔405) (檔406) (檔407) (檔408) (檔409) (檔410) (檔411) (檔412) (檔413) (檔414) (檔415) (檔416) (檔417) (檔418) (檔419) (檔420) (檔421) (檔422) (檔423) (檔424) (檔425) (檔426) (檔427) (檔428) (檔429) (檔430) (檔431) (檔432) (檔433) (檔434) (檔435) (檔436) (檔437) (檔438) (檔439) (檔440) (檔441) (檔442) (檔443) (檔444) (檔445) (檔446) (檔447) (檔448) (檔449) (檔450) (檔451) (檔452) (檔453) (檔454) (檔455) (檔456) (檔457) (檔458) (檔459) (檔460) (檔461) (檔462) (檔463) (檔464) (檔465) (檔466) (檔467) (檔468) (檔469) (檔470) (檔471) (檔472) (檔473) (檔474) (檔475) (檔476) (檔477) (檔478) (檔479) (檔480) (檔481) (檔482) (檔483) (檔484) (檔485) (檔486) (檔487) (檔488) (檔489) (檔490) (檔491) (檔492) (檔493) (檔494) (檔495) (檔496) (檔497) (檔498) (檔499) (檔500) (檔501) (檔502) (檔503) (檔504) (檔505) (檔506) (檔507) (檔508) (檔509) (檔510) (檔511) (檔512) (檔513) (檔514) (檔515) (檔516) (檔517) (檔518) (檔519) (檔520) (檔521) (檔522) (檔523) (檔524) (檔525) (檔526) (檔527) (檔528) (檔529) (檔530) (檔531) (檔532) (檔533) (檔534) (檔535) (檔536) (檔537) (檔538) (檔539) (檔540) (檔541) (檔542) (檔543) (檔544) (檔545) (檔546) (檔547) (檔548) (檔549) (檔550) (檔551) (檔552) (檔553) (檔554) (檔555) (檔556) (檔557) (檔558) (檔559) (檔560) (檔561) (檔562) (檔563) (檔564) (檔565) (檔566) (檔567) (檔568) (檔569) (檔570) (檔571) (檔572) (檔573) (檔574) (檔575) (檔576) (檔577) (檔578) (檔579) (檔580) (檔581) (檔582) (檔583) (檔584) (檔585) (檔586) (檔587) (檔588) (檔589) (檔590) (檔591) (檔592) (檔593) (檔594) (檔595) (檔596) (檔597) (檔598) (檔599) (檔600) (檔601) (檔602) (檔603) (檔604) (檔605) (檔606) (檔607) (檔608) (檔609) (檔610) (檔611) (檔612) (檔613) (檔614) (檔615) (檔616) (檔617) (檔618) (檔619) (檔620) (檔621) (檔622) (檔623) (檔624) (檔625) (檔626) (檔627) (檔628) (檔629) (檔630) (檔631) (檔632) (檔633) (檔634) (檔635) (檔636) (檔637) (檔638) (檔639) (檔640) (檔641) (檔642) (檔643) (檔644) (檔645) (檔646) (檔647) (檔648) (檔649) (檔650) (檔651) (檔652) (檔653) (檔654) (檔655) (檔656) (檔657) (檔658) (檔659) (檔660) (檔661) (檔662) (檔663) (檔664) (檔665) (檔666) (檔667) (檔668) (檔669) (檔670) (檔671) (檔672) (檔673) (檔674) (檔675) (檔676) (檔677) (檔678) (檔679) (檔680) (檔681) (檔682) (檔683) (檔684) (檔685) (檔686) (檔687) (檔688) (檔689) (檔690) (檔691) (檔692) (檔693) (檔694) (檔695) (檔696) (檔697) (檔698) (檔699) (檔700) (檔701) (檔702) (檔703) (檔704) (檔705) (檔706) (檔707) (檔708) (檔709) (檔710) (檔711) (檔712) (檔713) (檔714) (檔715) (檔716) (檔717) (檔718) (檔719) (檔720) (檔721) (檔722) (檔723) (檔724) (檔725) (檔726) (檔727) (檔728) (檔729) (檔730) (檔731) (檔732) (檔733) (檔734) (檔735) (檔736) (檔737) (檔738) (檔739) (檔740) (檔741) (檔742) (檔743) (檔744) (檔745) (檔746) (檔747) (檔748) (檔749) (檔750) (檔751) (檔752) (檔753) (檔754) (檔755) (檔756) (檔757) (檔758) (檔759) (檔760) (檔761) (檔762) (檔763) (檔764) (檔765) (檔766) (檔767) (檔768) (檔769) (檔770) (檔771) (檔772) (檔773) (檔774) (檔775) (檔776) (檔777) (檔778) (檔779) (檔780) (檔781) (檔782) (檔783) (檔784) (檔785) (檔786) (檔787) (檔788) (檔789) (檔790) (檔791) (檔792) (檔793) (檔794) (檔795) (檔796) (檔797) (檔798) (檔799) (檔800) (檔801) (檔802) (檔803) (檔804) (檔805) (檔806) (檔807) (檔808) (檔809) (檔810) (檔811) (檔812) (檔813) (檔814) (檔815) (檔816) (檔817) (檔818) (檔819) (檔820) (檔821) (檔822) (檔823) (檔824) (檔825) (檔826) (檔827) (檔828) (檔829) (檔830) (檔831) (檔832) (檔833) (檔834) (檔835) (檔836) (檔837) (檔838) (檔839) (檔840) (檔841) (檔842) (檔843) (檔844) (檔845) (檔846) (檔847) (檔848) (檔849) (檔850) (檔851) (檔852) (檔853) (檔854) (檔855) (檔856) (檔857) (檔858) (檔859) (檔860) (檔861) (檔862) (檔863) (檔864) (檔865) (檔866) (檔867) (檔868) (檔869) (檔870) (檔871) (檔872) (檔873) (檔874) (檔875) (檔876) (檔877) (檔878) (檔879) (檔880) (檔881) (檔882) (檔883) (檔884) (檔885) (檔886) (檔887) (檔888) (檔889) (檔890) (檔891) (檔892) (檔893) (檔894) (檔895) (檔896) (檔897) (檔898) (檔899) (檔900) (檔901) (檔902) (檔903) (檔904) (檔905) (檔906) (檔907) (檔908) (檔909) (檔910) (檔911) (檔912) (檔913) (檔914) (檔915) (檔916) (檔917) (檔918) (檔919) (檔920) (檔921) (檔922) (檔923) (檔924) (檔925) (檔926) (檔927) (檔928) (檔929) (檔930) (檔931) (檔932) (檔933) (檔934) (檔935) (檔936) (檔937) (檔938) (檔939) (檔940) (檔941) (檔942) (檔943) (檔944) (檔945) (檔946) (檔947) (檔948) (檔949) (檔950) (檔951) (檔952) (檔953) (檔954) (檔955) (檔956) (檔957) (檔958) (檔959) (檔960) (檔961) (檔962) (檔963) (檔964) (檔965) (檔966) (檔967) (檔968) (檔969) (檔970) (檔971) (檔972) (檔973) (檔974) (檔975) (檔976) (檔977) (檔978) (檔979) (檔980) (檔981) (檔982) (檔983) (檔984) (檔985) (檔986) (檔987) (檔988) (檔989) (檔990) (檔991) (檔992) (檔993) (檔994) (檔995) (檔996) (檔997) (檔998) (檔999) (檔1000) (檔1001) (檔1002) (檔1003) (檔1004) (檔1005) (檔1006) (檔1007) (檔1008) (檔1009) (檔1010) (檔1011) (檔1012) (檔1013) (檔1014) (檔1015) (檔1016) (檔1017) (檔1018) (檔1019) (檔1020) (檔1021) (檔1022) (檔1023) (檔1024) (檔1025) (檔1026) (檔1027) (檔1028) (檔1029) (檔1030) (檔1031) (檔1032) (檔1033) (檔1034) (檔1035) (檔1036) (檔1037) (檔1038) (檔1039) (檔1040) (檔1041) (檔1042) (檔1043) (檔1044) (檔1045) (檔1046) (檔1047) (檔1048) (檔1049) (檔1050) (檔1051) (檔1052) (檔1053) (檔1054) (檔1055) (檔1056) (檔1057) (檔1058) (檔1059) (檔1060) (檔1061) (檔1062) (檔1063) (檔1064) (檔1065) (檔1066) (檔1067) (檔1068) (檔1069) (檔1070) (檔1071) (檔1072) (檔1073) (檔1074) (檔1075) (檔1076) (檔1077) (檔1078) (檔1079) (檔1080) (檔1081) (檔1082) (檔1083) (檔1084) (檔1085) (檔1086) (檔1087) (檔1088) (檔1089) (檔1090) (檔1091) (檔1092) (檔1093) (檔1094) (檔1095) (檔1096) (檔1097) (檔1098) (檔1099) (檔1100) (檔1101) (檔1102) (檔1103) (檔1104) (檔1105) (檔1106) (檔1107) (檔1108) (檔1109) (檔1110) (檔1111) (檔1112) (檔1113) (檔1114) (檔1115) (檔1116) (檔1117) (檔1118) (檔1119) (檔1120) (檔1121) (檔1122) (檔1123) (檔1124) (檔1125) (檔1126) (檔1127) (檔1128) (檔1129) (檔1130) (檔1131) (檔1132) (檔1133) (檔1134) (檔1135) (檔1136) (檔1137) (檔1138) (檔1139) (檔1140) (檔1141) (檔1142) (檔1143) (檔1144) (檔1145) (檔1146) (檔1147) (檔1148) (檔1149) (檔1150) (檔1151) (檔1152) (檔1153) (檔1154) (檔1155) (檔1156) (檔1157) (檔1158) (檔1159) (檔1160) (檔1161) (檔1162) (檔1163) (檔1164) (檔1165) (檔1166) (檔1167) (檔1168) (檔1169) (檔1170) (檔1171) (檔1172) (檔1173) (檔1174) (檔1175) (檔1176) (檔1177) (檔1178) (檔1179) (檔1180) (檔1181) (檔1182) (檔1183) (檔1184) (檔1185) (檔1186) (檔1187) (檔1188) (檔1189) (檔1190) (檔1191) (檔1192) (檔1193) (檔1194) (檔1195) (檔1196) (檔1197) (檔1198) (檔1199) (檔1200) (檔1201) (檔1202) (檔1203) (檔1204) (檔1205) (檔1206) (檔1207) (檔1208) (檔1209) (檔1210) (檔1211) (檔1212) (檔1213) (檔1214) (檔1215) (檔1216) (檔1217) (檔1218) (檔1219) (檔1220) (檔1221) (檔1222) (檔1223) (檔1224) (檔1225) (檔1226) (檔1227) (檔1228) (檔1229) (檔1230) (檔1231) (檔1232) (檔1233) (檔1234) (檔1235) (檔1236) (檔1237) (檔1238) (檔1239) (檔1240) (檔1241) (檔1242) (檔1243) (檔1244) (檔1245) (檔1246) (檔1247) (檔1248) (檔1249) (檔1250) (檔1251) (檔1252) (檔1253) (檔1254) (檔1255) (檔1256) (檔1257) (檔1258) (檔1259) (檔1260) (檔1261) (檔1262) (檔1263) (檔1264) (檔1265) (檔1266) (檔1267) (檔1268) (檔1269) (檔1270) (檔1271) (檔1272) (檔1273) (檔1274) (檔1275) (檔1276) (檔1277) (檔1278) (