Dr.Rock Heat Therapy Pain Relief Gun

Product Introduction

Dr.Rock BioStone® Graphene Heat Therapy Pain Relief Gun, is an intelligent physical therapy device that combines traditional Chinese medicine (Dr. Rock BioStone®), Mobilization (heat therapy) and TuiNa (massage). Using the new graphene heat injection technology, after heating the stone, the far-infrared rays, negative ions and ultrasonic waves benefit to the human body, and the cells to produce the same frequency co-resonance, balancing the cellular potential, activation of the cells, promote blood circulation, relieve pain, relaxation of the muscles, and recovery of physical performance. It is effective in relieving sciatica, frozen shoulder, muscle aches and pains, and unblocking meridians throughout the body.

Six Key Features

- Heating 40~60°C
- Natural mineral BioStone
- Graphene heating
- Far infrared
- Five adjustment angles
- Retractable handle

Product Diagram

Instructions for use

1. Angle Adjustment

Press and hold the rotary button to rotate the following arm (all positions are available). Release the rotary button and allow the unit to settle. The arrow indicates that it is slightly clamped.

2. Length Adjustment

Unlock by pressing in the direction of the icon, applying force to extend or shorten the desired length, and then locking.

3. Heating Function

Switching on: Press and hold the power button and release when the power indicator appears. (Power on default: Start-up Heating Level 1 Temperature 40°C)
Switch off: Long press the power button and release it after the power indicator disappears.
Heating: Default Heating Level 1 at power on, short press the power button to cycle to Level 2 → Level 3 → off Heating.
In pause mode, the unit will switch off automatically after 3 minutes.

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Warranty, Cleaning, Storage

1. Before the product is used, clean it with a slightly damp towel and then dry it with a soft cloth.

2. Before and after using the BioStone-therapy head, wipe and disinfect the surface of the BioStone with medical alcohol.

3. Turn off the power after each use of the product.

4. When using the product for a long time or during travel, please place it in the provided special storage bag.

5. When the product is not in use for a long time, it should be charged every month to avoid battery damage or moisture.

6. When placing the product, pay attention to waterproofing, fireproofing, acid-proofing, and preventing impact or falling. Additionally, prevent mirrors from using it.

7. During long-term use of the product, it is important to maintain a good quiet and noise reduction experience, it is necessary to apply lubricant appropriately between the sleeve parts of the machine body (as shown in the diagram on the right).

Specifications

Model: H-MG300
Name: Dr.Rock BioStone® Graphene Heat Therapy Pain Relief Gun
Impact stroke: 12mm / Strength: 12KG
Gear adjustment: 1000 rpm ~ 3200 rpm
Heat treatment: 1 gear 40°C, 2 gear 50°C, 3 gear 60°C
Power input: USB-C 2A/5V
Battery: 1600mAh, 14.8V lithium battery
Operating Temperature Range: 0°C ~40°C
Endurance Time: 2-8 hours / Charging time: about 5.5 hours
Product size: 98*239*81mm (Folded) 148*330*81mm (90° extension)
Main unit weight: 690g / Product color: Black

Package content

- Dr.Rock Heat Therapy Pain Relief Gun x1
- USB-C Cable x1
- BioStone tapered head x1
- Storage Bag x1
- User's manual x1

BioStone Far-Infrared Heat Therapy Principle

Warming mobilization therapy is an important method in traditional Chinese medicine (BioStone therapy). It originated from placing heated Chinese stones on affected areas or using them for massage, which can alleviate pain and treat ailments. This therapy is known as the warming mobilization of BioStones.

Electrothermal BioStone therapy is a kind of warm mobilization therapy, i.e., the electric current will be heated through a special instrument to BioStone, the far-infrared wave radiation heat produced by BioStone heat is transmitted in the meridians to achieve the purpose of treating various diseases, is the clever combination of traditional BioStone therapy and modern high-tech. At the same time, BioStone has a storage energy field, the role of the human skin surface can produce far infrared radiation, its bandwidth is extremely wide, far infrared frequency up to the range of 7-20 microns. BioStone every friction of the human body can produce beneficial to the health of the ultrasonic pulse, the average number of ultrasonic pulses up to 3708 times, the frequency range of 2-2 million Hz is the highest of all materials.

Massage Impact Therapy

Impact on Skin and Subcutaneous Tissue:
Massage impact therapy influences the skin and underlying tissues by loosening thickened connective tissue and fascia. It breaks down barriers that hinder the release of fascia contractions, removes adhesions and internal scar tissue after injury or surgery, and improves venous return and lymphatic circulation. It also enhances muscle extensibility, especially in cases of muscle spasms, stiffness, and limitations. The therapy stimulates reflexive contractions of muscle fibers and increases the range of motion by relieving joint stiffness.

Impact on the Central and Sympathetic Nervous System:
Massage impact therapy stimulates local muscle receptors to alleviate imbalance between the parasympathetic and sympathetic nervous systems, and promotes relaxation.

Impact on Muscle, Fascia, and Joint Pain:
The most common causes of muscle pain are trauma, excessive tension, or tightness. After injury or trauma, impact therapy accelerates the healing process in inflamed areas. Similarly, when certain parts of the body are overused or tense, impact therapy enhances blood circulation to the affected areas and alleviates potential tension.

Massage impact therapy aids in recovering the body from the stress of intense exercise. It helps regulate and rebuild training by eliminating muscle fatigue, nourishing cells, restoring tissue elasticity, and maximizing muscle relaxation. The therapy stimulates trigger points and accelerates wound healing.

Recovery after Surgery and Injury

Following trauma of surgical repair, impact therapy provides stimulation during the body's inflammatory healing process. It helps prevent the formation of scar tissue and adhesions.

Relaxation:
Massage impact therapy relaxes the body, promotes comfort, and enhances overall well-being, thereby increasing the patient's satisfaction.

Indications and Contraindications

The Dr.Rock series of deep penetration therapy devices can effectively alleviate pain, bursitis, muscle spasms, neuritis, and other musculoskeletal conditions caused by joint sprains, overexertion, and abnormal posture. BioStone can be used to treat the scar tissue involved in most of these conditions. When using Dr.Rock, it is important to avoid applying excessive pressure and engaging in prolonged therapy sessions. Additionally, sensitive areas such as the head, face, and cervical vertebrae should be avoided. If there is inflammation, swelling, or an increase in pain, therapy should be immediately discontinued. Current research is focused on studying the effects of vibration training on specific medical conditions. This research has the potential to reduce the number of contraindications mentioned earlier. Based on practical experience, combining vibration training with a physiotherapy program is often recommended in many cases. However, it is important to seek advice from a doctor, expert, or physiotherapist and perform the therapy under their guidance or supervision.

Indications

- ✓ Pain and spasms caused by muscle contusions, sprains, and strains.
- ✓ Assisting in the flow of edema fluid in swollen areas.
- ✓ Relieving thickened connective tissue and fascia.
- ✓ Reducing lactic acid accumulation in muscles.
- ✓ Increasing joint mobility.
- ✓ Increasing muscle fatigue.
- ✓ Sciatica (BioStone hot compress physiotherapy).
- ✓ Sacrospinal periarthritis (BioStone hot compress physiotherapy).
- ✓ Patellaritis (BioStone hot compress physiotherapy).

Contraindications

- ✗ Aneurysm, hemorrhage, and use of blood thinners.
- ✗ Heart disease with a cardiac pacemaker or defibrillator.
- ✗ Pregnancy and cancer.
- ✗ It should not be used within 3 inches of internal fixation when the joint internal fixation is less than 90 days.
- ✗ Sensitive areas: head, face, cervical spine, spinal vertebrae.
- ✗ Locations with less muscle coverage near bones, bony ball, foot, and back of the hand.

Abnormal Malfunctions and Troubleshooting

Issue Description	Possible Cause	Recommended Actions
The product is not turning on.	1. Insufficient battery power. 2. The product has been operating continuously for 15 minutes. 3. The protection mechanism is triggered.	Change the product (if not used for a long time, charge for more than 7.5 hours). Restart the product as it may have entered long-term use protection mode. Restart the product.
The product automatically shuts down during use.	1. Low battery. 2. The product has been operating continuously for 15 minutes. 3. The protection mechanism is triggered.	1. Charge the product. 2. Restart the product as it may have entered long-term use protection mode. 3. Restart the product.
The product cannot be charged.	1. Check if the port is properly connected and ensure that the adapter supports 5V/2A. 2. Changing is not possible if the product temperature or ambient temperature exceeds 40°C/104°F.	1. Try to replace the adapter or plug it into a compatible one. 2. If the product temperature is high, wait for it to cool down or recharge it in a normal temperature environment.
The therapy head falls off during usage.	Check if the therapy head is installed properly.	Refer to the operation reference manual for instructions on how to install the therapy head.
The temperature of the BioStone head is too high.	Turn off the product and wait for the temperature of the BioStone to cool down. Alternatively, you can lower the intensity of the hot compress mode before use.	Contact after-sales.
Excessive noise during product use.	The partial slits may have dry grease.	Contact after-sales.

Physiotherapy Considerations

After Injury

The first physiotherapy session after injury should gradually progress based on the patient's condition and comfort.

- Range of motion (ROM) assessment.
- Determine the muscles and areas to target based on the physiotherapy plan.
- Perform gentle and systemic physiotherapy.
- Apply pressure equivalent to the weight of the device itself.

Lymphatic Drainage

- Perform gentle and systemic physiotherapy.
- Move slowly from the distal end to the proximal end of the limb, following the direction of venous return to the heart.

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